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
The Stefan-Boltzmann Lamp is equipped with a filament of tungsten.

The resistance of the filament varies with temperature. Therefore we can determine the wire temperature while the lamp is in operation by measuring this resistance. The resistance is determined using Ohm's Law.

To avoid measuring the resistance in wires, plugs, sockets, etc., it is important to measure the voltage as close to the filament as possible. The lamp is equipped with a separate set of connectors (the upper pair) to allow a voltmeter to be connected directly to the lamp socket.

Emission of light: Stefan-Boltzmann's law

The Stefan-Boltzmann Law states that the power P_R radiated from a black body with surface area A and absolute temperature T is given by

$$P_R = A \cdot \sigma \cdot T^4$$

Propagation of Light: The inverse-square law

The radiation intensity (power per area) I_R from a point source decreases with distance r like

$$I_R = \frac{P_S}{4 \cdot \pi \cdot r^2}$$

Temperature dependency of tungsten

From the actual resistance R of the filament and R_{300} (the resistance at 300 K) its *absolute* temperature T is found by first calculating the relative resistance change $x = R/R_{300}$ which is then inserted into a sixth order polynomial:

$$T = a_6 x^6 + a_5 x^5 + a_4 x^4 + a_3 x^3 + a_2 x^2 + a_1 x + a_0$$

The coefficients are given by:

a_6 :	$-9.915 \cdot 10^{-5}$	a_2 :	-19.76
a_5 :	$6.794 \cdot 10^{-3}$	a_1 :	267.3
a_4 :	-0.1840	a_0 :	52.38
a_3 :	2.520		

The following formula gives us R_{300} from a measurement performed at room temperature t_{REF} (measured in °C)

$$R_{300} = R_{\text{REF}} \cdot \left(1 + \frac{26.8^\circ\text{C} - t_{\text{REF}}}{208.3^\circ\text{C} + t_{\text{REF}}} \right)$$

From our web site you can download a spreadsheet with these formulae (search for 277022).

Experiments

You can download complete lab manuals from our web site www.frederiksen.eu:

133720 Stefan-Boltzmann's law

133710 The inverse-square law

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© Frederiksen Scientific A/S