Bredbåndssensor

12.04.13

2872.81 AC

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2872.81 AC

Description

This detector is sensitive to electromagnetic radiation in the wavelength region from 250 nm to 20 μm , i.e. from the ultraviolet through the visible and near infrared to the far infrared region of the spectrum. The far infrared from 3-20 μm included heat radiation e.g. from the human body.

Although the sensitivity is less than for silicon detectors, this detector type has the important advantage of having a very flat spectral responsivity curve. The sensitivity is practically wavelength independent, depending only on the incident power per unit area. The sensor is calibrated (see below).

Necessary accessories

The sensor may be connected directly to e.g. Pasco's PS-2158 analog adaptor or a Pasco 500 interface.


Another option is to connect the sensor to a 2515.60 battery box.

Mounting, operation and maintenance

The sensor is connected to the read out unit and then directed towards an object which emits electromagnetic radiation in the 250 nm - 20 μm region. The response time is about 0.5 s (to 1 %) Avoid touching the sensitive detector window, and avoid moisture. As the sensor window consists of KBr which is hygroscopic, the sensor must be stored in a closed plastic bag together with some silica gel.

Specifications

Dimensions:	$\varnothing 25 \times 100 \text{ mm}$
Sensor element:	2 x 2 mm
Sensitive angle:	38° (full sensitivity – see figure)
Output signal:	0 to 5 V
Precision:	3 %

Calibration

A calibration factor is specified on a label on the sensor. The voltage (unit: V) is multiplied with this factor, giving the intensity (unit: W/m²)

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© A/S Søren Frederiksen, Ølgod