Instruction Manual

1. PRODUCT DESCRIPTION

This pipette is a general purpose pipette for sampling and dispensing accurate liquid volumes. The pipettes operate on the air displacement principle and use disposable tips. The black tip cone can be sterilized (121°C,0.15Mpa).

The pipettes cover a volume range of 0.5ul to 5ml

All pipettes have been quality tested according to our company's demands. The quality control according to our company's involves gravimetric testing of each pipette with distilled water at 22°C using the manufacter's original tips.

1.1. Fixed volume pipette

Cat. No.	Volume(ul)	Inacc,%	Impr.%
A3001	0,1:	±4,8	≤13.01
A3002	0.2	±4,6	≤12.81
A3003	0.25	±4,6	≤ 2,8
A3004	0.5	±4.5	≤ 2.81
A3005	11	±4.2	≤ 2,8
A3006	2	±4,0	≤12,61
A3007	2,5	±4.0	≤12,61
A3008	5	±3_8	≤12.41
A3009	10	± 3_0	€12_01
A3010	20	±2,4	≤ 1.01
A3011	25	±2,4	€11.01
A3012	50	± 2, 4	≤11.01
A3013	100	±1,5	≤10.81
A3014	200	±1;2	≤ 0_7
A3015	250	±1,2	≤10.71
A3016	500	±1,2	≤+0.71
A3017	1000	±1,0	≤10.51
A3018	2000	± 1,0	≤10.51
A3019	2500	±1,0	≤10.51
A3020	5000	±1.0	≤10.51

1.2.Adjustable volume pipette

Cat No	Voleme Rangeint)	Increment (ut)	Test Volunie (ui)	Inacc (%)	Impr (%)	Tip ID
А3021 0.1-2.5µl	0.1-2.5µl	0. 05 μ!	2.5µI	±3.0	≤ 3.0	10µ1
			1.25pl	±6.0	≤ 14.0 [
		0,25µl	±12.0	≤ 6,0		
A3022 0,5-10µ			10 ul	±3.0	≤ 12,01	
	0.5-10µl	0.5-10µl 0.1µl	5 Д	_ ±6.0	≤ 13,01	10 01
			14	±10.0	≤ 3,01	
A3023 2-20µl			20 M	±3.0	≤ 1131	
	2-20µl	0.5 µl	10 pl	±3.0	≤ 2,01	200µ1
		2 µ]	±10	≤ 3,01		
A3024	5-50µI 0,5 µI		50 µl	±2-4	≤ 11.51	200µl
		0,5 µl	25 ⊭l	±2.8	≤ 2,0	
			3st	±6.0	≤ 13,01	
			100 µl	±l.ā	≤ 11.01	

Cat No	Volume Range(m)	Increment	Yest Volume	Inace (%)	1mpt (%)	Tip ID
A3026 20-200µl		200 #	±1.5	≤ 11.01		
	20-200ul	1 µl	100 14	±13	≤ 11.01	200µl
	20 200 51		20 H	±3.0	≤ [1.5]	
			200 pl	±13	≤ 11.01	200µ1
Λ3027 50-200μl	50-200ul	1 pi	100 (4	±1.5	≤ 1.01	
	30 200		50 M	±22	≤ 11.51	
			1000 pl	±1.2	€ [0.5]	1000 µ
А3028 100-1000µ	100-1000ul	5 µl	500 A	±1.5	≤ 10.51	
	100 1000		100 pi	±1.5	≤ 1,01	
		5 µl	1000 pl	±1.2	≤ 10.51	1000 pl
A3029 200-10	200-1000µl		500 µl	±1.5	≤ 10.51	
			200 회	±1.5	< 11.01	
A3030 1000-5000µ		50 ul	5000 4	±1.0	≤ [0.5]	5000 pl
	1000-5000ul		2500 H	±1,0	≤ 10.51	
			1000 st	±1.2	< 10.51	
			10ml	±0.5	≤ 10.21	10ml
A3031	2-10ml	0.1 ml	5ml	±0.8	≤ 03	
A3031 2 101			2m1	±1.0	≤ 10.31	
A3045 8CH0.5-10			10µl	±3_0	< 12.01	10μ1
	8CH0.5-10µl	0 l µl	5µl	±6.0	≤ 3.01	
			1 μ1	±10.0	≤ 13.01	
		ι 0,5 μ1	50µl	±2,4	€ 11.51	200µl
A3046 8C115-50	8C115-50µl		25µl	±2.8	€ 12.01	300µl
			5µl	±6.0	≤ 3.0	
			300µ1	±1.8	≤ 11.51	300µl
A3047	8CH50-300µl	5 µi	150µl	±2.0	≤ 2.01	
AJU47 OCID	oc.iii		50µ1	±2.2	≤ 13.01	
A3061 12C110.5-1			10µ1	±3.0	≤ 2.0	10µ1
	120110.5-10	0.1 pl	5µl	±6.0	≤ 13.01	
	124.1015.101	attorio a color	1µl	±10.0	€ 13.01	
A3062 12		0,5 μ1	50µl	±2.4	€ 11.51	200μ 300μ
	12CH5-58pl		25µl	+2.8	≤ 12.01	
	120.02-200		5µl	±6.0	≤ 13.01	
_	-	5 ді	300ul	±1.8	€ 11.51	300μ
A3063	12C1150-300µ		150µ1	±2.0	≤ 12.01	
A3063			50µI	±2.2	≤ 3.0	

1.3.Tips

These detachable, disposable tips are made of natural colour polypropylene. Tips packed in bulk or box

Note: never pipette liquid without attaching a tip to the pipette.

1.4. Black tip cone

The black tip cone can be sterilized(121°C,0-15Mpa)



Fig.1

1.5. Pipette materials



2. PACKAGE

The pipette package contains the following items:

- ♦ Pipette
- ♦ Calibration/opening tool
- ♠ Grease
- ♦ Instructions for use

3.PIPETTE OPERATION

3.1. Volume setting

The volume of the pipette is clearly shown through the handle grip window. The delivery volume(variable volume pipettes only) is set by turning the thumb button clockwise or anticlockwise(Fig.3). When setting the volume, please make sure that:



- ♦ the digits are completely visible in the display window.
- the selected volume is within the pipette's specified range.

Note: Using excessive force to turn the push button outside the range may jam the mechanism and damage the pipette.

3.2. Sealing and ejecting tips

Before fitting a tip, make sure that the pipette tip cone is clean. Press the tip on the cone of the pipette firmly to ensure an airtight seal. The seal is tight when a visiblesealing ring forms between the tip and the black tip cone(Fig.4).

Each pipette is fitted with a tip ejector to help eliminate the safety hazards associated with contamination. The tip ejector needs to be pressed firmly downwards to ensure proper tip ejection(Fig.5). Make sure that the tip is disposed of into a suitable waste container.

3.3.Removing the tip ejector

Please removing the tip ejector collar by following instructions (Fig.6)

- ♦ Push the tip ejector down: (1)
- Push the opening tool pin between the ejector bar and ejector collar to release the locking mechanism. (2)
- ♦ Pull the ejector collar off (3)
- Pushing the collar locking pin into the ejector bar hole, keeping the tip ejector pushed down.



Fig.3



Fig.4



Fig.5



4.PIPETTING TECHNIQUES

4.1 Forward technique

For the best possible accuracy, please push and release the button slowly at all times, particularly when working with liquids that have a high viscosity, Never allow the operating button to snap back. Ensure that a clean tip is firmly pushed on to the tip cone of the pipette and there are no foreign particles in the tip itself.

Hold the pipette vertically during aspiration, Make sure that the liquid and container vessel are clean and that the pipette, tips and the liquid are at the same temperature,

1. Depress the operating button to the first stop(Fig.8B).

2. Dip the tip under the surface of the liquid about 2-3mm deep and slowly release the operating button (Fig. 8A). Withdraw the tip from the liquid touching it against the edge of the container to remove excess liquid.

3. Deliver the liquid by gently depressing the operating button to the first stop (Fig. 8B). After a delay of about one second, continue to depress the operating button all the way down to the second stop (Fig. 8C). This action will empty the tip.

4. Release the operating button to the ready position. (Fig. 8A). If necessary, change the tip and continue with the pipetting.



Starting position Fig.8A



First stop Fig.8B



Second stop Fig.SC

4.2 Reverse technique

The reverse technique is suitable for dispensing liquids having high viscosity or a tendency to foam easily. The technique is also recommended for dispensing very small volumes. This is achieved by filling and emptying the tip.

1. Depress the operating button all the way down to the second stop (Fig. 8C).
2. Dip the tip under the surface of the liquid about 2-3mm deep and slowly

release the operating button (Fig. 8A). This action will fill the tip. Withdraw the tip from the liquid touching it against the edge of the container to remove excess limid.

3. Deliver the preset volume by gently depressing the operating button to the first stop (Fig.8B). Hold the operating button at the first stop. Some liquid will remain in the tip and should not be included in the delivery.

4. The remaining liquid is either discarded with the tip or pipette back into the container.

- Prewet tip before aspirating the liquid by filling and emptying the tip 5 times. This is important especially when dispensing liquids which have a viscosity and density different from water.
- Always control the push button movements with the thumb to ensure consistency.
- When pipetting liquids at a temperature different from ambient, prewet the tip several times before use.

6. STORAGE

When the pipette is not in use make sure it is safely stored in a vertical position.

7. CALIBRATION

Each pipette has been factory-tested and certified at 22°C according to ISO 8655/DIN 12650. The Fmax user should not exceed the Fmax by more than 100%. Note: Pipette specifications are guaranteed only with manufacturer's tips.

7.1 Checking calibration

Weighing should take place at 20-25°C, constant to 0.5°C. Avoid drafts

- Set the desired testing volume of your pipette.
- Carefully fit tip onto the tip cone.
- Prewet tip with distilled water by pipetting the selected volume 5 times.
- Carefully aspirate the liquid, keeping the pipette vertical.
- Pipette distilled water into a tared container and read the weight in mgs. Repeat at least five times and record each result. Use an analytical balance with a readability of 0.01 mgs. To calculate the volume, divide the weight of the water by its density(at 20°C: 0.9982).
- Calculate the F-value by using the following equation;
- ◆ F= | inaccuracy(ul) | +2×imprecision(ul)
- Compare the calculated F-value to the corresponding Fmax user. Otherwise check both your accuracy and precision, when necessary, proceed to recalibration procedure.

7.2 Recalibration

I.Place the calibration tool into the holes of the calibration adjustment lock(under the operating button)(Fig. 9).

2.Turn the adjustment lock anticlockwise to decrease and clockwise to increase the volume.
 3.Reapeat performance test procedure form step I until the pipetting results are correct.



To maintain the best results from your pipette each unit should be checked every day for cleanliness. Particular attention should be paid to the tip cones.

The pipettes have been designed for easy in-house service. However, we also provide complete repair and calibration service. Please return your pipette to your local representative for repair or recalibration. Before returning, please make sure that it is free from all contamination. Please advice our Service Representative of any hazardous materials which may have been used with your pipette.

Note: Check the performance of your pipette regularly e.g.every 3 months and always after in-house service or regularly.

8.1. Cleaning your pipette

To clean your pipette, use ethanol and soft cloth or lint-free tissue. It is recommended to clean the lip cone regularly.

8.2.In-house maintenance

1. Hold down the tip pipette.

2.Place the tooth of the opening tool between the tip ejector and the tip ejector collar to release the locking mechanism. (Fig.10)

3.Carefully release the tip ejector and remove the ejector collar.
4.Place the wrench end of the opening tool over the tip cone, turning it

anticlockwise. Do not use any other tools . (Fig.11)

5. Wipe the piston, the O-ring and the tip cone with ethanol and a lint-free cloth.

Note: Models up to $10\mu l$ have a fixed O-ring locked inside the tip cone.

Therefore, the O-ring cannot be accessed for maintenance.

6.Before replacing tip cone it is recommended to grease the piston slightly using the silicone grease provided.

Note: Excessive use of grease may jam the piston.

7. After reassembling use the pipette (without liquid) several times to make

sure that the silicone grease is spread evenly. 8.Check the pipette calibration.





9. TROUBLE SHOOTING

The accompanying table is a guide to possible problems and their solutions.

Trouble	Possible cause	Solution	
Droplets left inside the tip Leakage or pipetted volume too small	Unsuitable tip Non-uniform wetting of the plastic Trip incorrectly attached Unsuitable tip Foreign particles	Use original tips Attach new tip Attach firmly Use original tips Clean the tip cone, attach new tip Clean and grease O-ring and	
	Instrument contaminated Insufficient amount of grease on piston and O-ting O-ting not correctly positioned or damaged Incorrect operation	piston, clean the tip cone Grease accordingly Change the O-ring Follow instructions carefully Recalibrate according to	
Push button jammed or moves erratically	Piston contaminated Penetration of solvent vapours	Clean and grease O-ring and piston, clean the tip cone	
Pipette blocked, aspirated volume too small	Liquid has penetrated tip cone and dried	Clean and grease O-ring and piston, clean the tip cone	
Tip ejector jammed or moves erratically	Tip cone and/or ejector collar contaminated	Clean the tip cone and the ejector collar	

10.WARRANTY INFORMATION

The pipettes are warranted for one year against defects in materials and workmanship. Should it fail to function in any period of time, please contact your local representative immediately. The warranty will not cover defects caused by normal wear or by using the pipette against the instructions given in this