Apparatet består af følgende dele:

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
11.02.15

1705.00 AC


The apparatus comprises the following:

- A generator
- The turbine wheel made of black plastic
- The turbine housing made of transparent acrylic
- A cover plate

Introduction

The apparatus comprises a generator mounted with a turbine wheel with 8 cups. A strong water or air jet will cause the turbine wheel to rotate and make the generator produce electricity. To avoid splashes of water, cover the turbine housing with the cover plate and apply the water through the hole. Connect a load (e.g. a light bulb) to the generator terminals, which are available as 4 mm security sockets.

Each cup is shaped to achieve the highest effect when a high-pressure water jet hits the cup in one end to be returned in the other (see diagram). In this way, the highest momentum change of the water jet is achieved. In theory the cup will get the same change in momentum but in opposite direction. Therefore, a narrow high-pressure tap is advised e.g. mount "Water jet vacuum pump, PP" (Item no.: 0690.20). Under these circumstances the turbine can deliver 2.4 V and 0.6 A which means that a bulb 2,5 V 0,3 A can be used (item no. 4250.10).

If an ordinary tap is used, you cannot expect more than a glow with a 1.5 volt, 90 mA bulb (item no.: 4250.05) - even at the highest water capacity. Therefore, we recommend using the turbine with the "Water jet vacuum pump, PP".

Technical specifications

Voltage max of generator:	12 V DC
Power:	up to 1.5 W (with "Water jet vacuum pump, PP")
Turbine Wheel:	140 mm diameter
Dimensions:	195 x 170 x 100 mm
Weight:	610 g

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© A/S Søren Frederiksen, Ølgod