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

Instruction manual for screen for sidemount. no. 1955.80


13.12.10

AC 1955.80

$s = 20 \text{ mm}$

Total width = 30 mm

Trigger pulse from the photocell (trigging on the rising edge)

The screen no. 1955.80 is intended for measuring the velocity of the gliders moving on the air track no. 1950.10 (with C-profile).

The screen mounts on the side of the glider by pressing it down over the steel pin protruding from the glider. Besides the steel pin a horizontal guiding rib matching the notch in the screen is used to fix the screen in a stable position. To avoid asymmetric load of the glider a mass of 5 g is supplied with the screen. Place the extra mass on the steel pin protruding from the gliders opposite side.

When measuring the gliders velocity you measure the passage time for a fixed distance e.g. the width of the screen. Start respectively stop of the time measuring unit (e.g. Scaler-Timer no. 2001.00 or 2001.50) may be done by means of a photocell unit (e.g. no. 1975.10).

When using more than 1 photocell unit for these measurements, even small variations in the photocell units light intensity may lead to faulty results, as the screen width will not be perceived identical by the photocells.

The screen no. 1955.80 has been designed especially with this problem in mind. The screen has an incision which means that the screen actually comprises of 2 screens instead of one. This design is able to eliminate the source of fault described above. When connected to a Scaler-Timer e.g. no. 2001.00 or 2001.50 measurement of time will start when the photocells light beam is interrupted by the first edge of the screen (rising flange) and stop when the second edge of the screen interrupts the light beam. I.e. the Timer-Scaler will start respectively stop time measurements when the photocell changes from light to dark. The screen is symmetric.