

2018-10-16 / HS

Produktmanual AC198510

English
manual
overleaf

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

Product manual AC198510

Dansk
manual:
Vend!

The switch box is used to send a start signal to a counter, an electronic stopwatch, or a SpeedGate when the current through a coil is interrupted. This will typically be useful during acceleration experiments on an air track. For that purpose a 195210 (or 195200) Electric launcher is used.

Application

Connect an 8 V DC power supply to the sockets labelled **+8 VDC** and **0** (See below.) Be careful not to reverse polarity. The electronics in the switch box may under rare circumstances be damaged by a voltage lower than 8 V – if you need a lower current, use for example a coil with a higher resistance (more windings).

Connect the coil from the launcher to the sockets marked **Coil**. The yellow **Hold** LED is lit when a voltage is applied across the coil. In this setting, the current in the coil will typically be higher than the coil's maximum rating, which is acceptable for short time intervals. When not needed, the current in the coil should be turned off (LED off). The **Coil** output must not be short circuited.

An electronic stopwatch or a universal counter (like 200280 or 200250) is connected to the sockets marked **Timer Signal**. A SpeedGate must be connected through a 198512 Signal limiter. See pictures above.

With a voltage across the coil, the signal to the stopwatch is *low* (approximately 0 V). The moment the current through the coil is interrupted, the signal to the stopwatch goes *high* (8-12 V). A signal limiter will reduce the voltage swing to 0-5 V.

The power supply

Voltages outside of parentheses apply to the new electric launcher 195210. Voltages in (parentheses) apply to the old model 195200.

With an unstabilized power supply without display (like 361050 or 361060), simply set the button at 12 V (8 V).

With an unstabilized power supply with display (like 361860), adjust the voltage to 12 V (8 V) with current running through the coil. The voltage will rise a little when current is shut off.

With a stabilized power supply (like 364000, 361600, 161700), adjust the voltage to 12 V (8 V).

The power supply must be able to provide the necessary current to the coil without activating any thermal circuit breaker or current limiter.

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© Frederiksen Scientific A/S