Beskrivelse

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
Vertical Acceleration Demonstrator 1990.00

2015-01-09

Ac 1990.00

Description

This device demonstrates the independence of the vertical and horizontal components of motion.

The demonstrator sets two steel balls in motion simultaneously. One ball falls vertically with no velocity in the horizontal direction. The other ball falls vertically but also has an initial horizontal velocity. It follows a parabolic trajectory. Your students will hear simultaneous “clicks” as the balls hit the floor at the same time.

The steel balls are stored in the ball storage compartment. Rotate the cover to remove the balls.

Procedure

1. Place the Vertical Acceleration Demonstrator on a level table, with the launch arm of the demonstrator beyond the edge of the table as shown. Be sure that the floor is also level.
2. Pull the launch arm back until it catches on the trigger.
3. Mount a steel ball in the small hole at each end of the arm.
4. Press the Trigger Knob down. One ball will drop straight down and the other will be projected horizontally.
5. Listen for the click of the balls as they strike the floor. If the apparatus and the floor are both level, the two balls will strike the floor simultaneously.

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© A/S Søren Frederiksen, Ølgod