
2017-06-13 / HS

Produktmanual AC199560 p. 1/4

English  
manual  
on p. 3

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
Product manual AC199560 p. 3/4

Dansk  
manual:  
Se p. 1

The ball launcher is used for studying projectile trajectories.

The equipment comprises a gun barrel with a clinometer, a mounting bracket, a ramrod, and three plastic balls.

#### Necessary equipment

The bracket is fastened to a table edge with a regular screw clamp.

The impact is recorded on a piece of paper using carbon paper.

We can recommend the following:

001510 Screw clamp

199215 Carbon paper A4, 100 sheets

#### Optional accessories

The ball launcher can be used with photogates from Pasco, Vernier or Fourier.

We have developed a bracket for mounting one or two photogates in front of the launcher.

199570 Photogate bracket

## Using the ball launcher

The launcher is typically used with a large (i.e. long) table. Two smaller tables of the same height can also be used. The horizontal plate of the bracket is fixed to the edge of the table with a screw clamp so that the cannon can be placed outside the table top with the muzzle towards the table.

The innermost hole in the bracket is only used for mounting the gun barrel if you don't use a photogate bracket.

The next hole makes room for the photogate bracket, even at firing angles close to horizontal. (In the photo, the barrel is mounted in this hole.)

The outermost hole can be used for horizontal firing of the ball away from the table. The gun barrel is then moved on the opposite side of the bracket. The screw clamp should possibly be angled slightly to make room.

The elevation of the gun barrel can be read on the integrated clinometer.

For all three mounting modes, *the bottom of the ball is at the height of the table at launch* - regardless of the firing angle.

The cannon is loaded with a plastic ball, by pressing it down the barrel with the ramrod. The blue trigger arm will engage the spring mechanism at three specific locations, corresponding to three different muzzle speeds.


Never look down the cannon to see if the launcher is loaded!

The launcher is fired by pulling the trigger arm using the cord.

After a test launch, you know approximately where the ball will land. Place a piece of A4 paper there. It can with advantage be taped to the table. A piece of carbon paper is placed on top of it. Then the firing is repeated and the bullet will now leave a mark on the paper at the impact point.

## Photogate bracket

This bracket is an optional accessory. The bracket enables you to place a photogate immediately beyond the muzzle of the launcher and another one 10 cm further away.

The position of the photogates can be adjusted lengthwise; the distance between them is fixed.

## Specifications

At a launch angle of 45°, horizontal ranges of approximately 1.1 m, 2.3 m and 4.5 m, respectively, are achieved for the three different muzzle speeds. With this launch angle, the ball will reach the following approximate maximum heights: 0.3 m, 0.6 m and 1.2 m. The corresponding maximum heights for vertical firing are approx. 0.6 m, 1.2 m and 2.3 m.

Note that these parameters are neither checked nor guaranteed; they are only specified here to give an impression of the space requirements of the setup.

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© Frederiksen Scientific A/S