
2018-10-03 / HS

Produktmanual AC200510 p. 1/4

English  
manual  
on p. 3

#### Timer til papirstrimmel.

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
### Ticker tape timer

Connect the timer to an AC voltage via the safety sockets. The current passes first a diode, then a magnet coil.

If you live a place with 50 Hz mains frequency, the diode lets current pass through the coil 50 times per second.

(And if your AC frequency is 60 Hz, the coil is energized 60 times per second.)

For each current pulse, a pin strikes a piece of carbon paper, making a dot at the paper tape. The distance between the dots therefore represents a time interval of 0.02 s ( 0.0167 s for 60 Hz ).

The timer is provided with a switch that must be pressed for the timer to work.

### Equipment needed

Ticker timer carbon paper, 35 mm, 50 pieces (NL121256)

Ticker tape, 18 mm (200520)

Power supply, capable of 12 V/2 A AC (for instance 361055, 361600 or similar)

Lab leads (e.g. 105753 – 2 pcs.)

For free fall experiments, one of these weights is used:


200550 Timer weight, 0.5 kg

200560 Timer weight, 0.25 kg

### Preparing the tape and the carbon paper

The ticker tape must pass through the two square hoops at the sides of the timer, and must likewise pass under the carbon paper.

New carbon paper is mounted by pushing it in beneath the clear plastic guide and coax it onto the steel pin.

Change the carbon paper when it is so worm that the dots are no longer prominent or if it starts to wrinkle.

### Adjustments

The armature with the pin for the dots sits in a flexible suspension. If needed, the spring can be tightened or loosened. The two finger nuts must be tightened against each other.

## Acceleration (due to gravity)

To measure the acceleration due to gravity, place the timer high on a stand rod that is held by a table clamp. Protect the floor against the impact of the weight. Use lab leads that are long enough that they can be out of the way.


Cut a length of ticker tape that is almost as long as the distance from the timer to the floor. Pass the tape through the timer and attach a weight at the end.

Hold the tape tight and bring the weight to rest.

Start the timer by pressing the switch – release the tape, so that the weight pulls it all the way through the timer – release the switch again.

Note: Other kinds of constantly accelerated motion can be investigated as well. You can e.g. fix the end of the tape to a cart that runs down an incline.

The ticker tape will look like this:


As we assume that the acceleration is constant, we don't need to decipher the very first millimeter of the tape, but can just define  $t = 0$  s in a clear dot close to the starting point.

### 50 Hz procedure

Mark every *fifth* dot, and calculate the speed of the weight as

$$v = \frac{\Delta s}{\Delta t}$$

$\Delta s$  is found as shown in the drawing here:

$\Delta t$  is  $2 \cdot 0.02$  s = 0.04 s.

### 60 Hz procedure

Use instead every *sixth* dot. The time for the *sixth* dot will then be 0.1 s as in the drawing above, but the time interval  $\Delta t$  is different:

$\Delta t$  is  $2 \cdot 0.01667$  s = 0.03333 s.

By plotting corresponding values of  $t$  and  $v$  in a coordinate system,  $g$  can be found as the slope of the best straight line through the data points.

## Conservation of mechanical energy

The tape let us not only measure the speed  $v$ , but also how far the weight has fallen  $L$ . This will let us find both the kinetic and the potential energy of the weight.

The zero point for potential energy can be chosen freely. If we define the potential energy of the weight to be zero when it is placed on the floor, the potential energy at the start point (at height  $h$ ) is given by

$$E_0 = m \cdot g \cdot h$$

The weight starts with kinetic energy zero.

The two energies are therefore found as

$$E_{\text{POT}} = E_0 - m \cdot g \cdot L$$

$$E_{\text{KIN}} = \frac{1}{2} \cdot m \cdot v^2$$


The sum of these is called the mechanical energy

$$E_{\text{MECH}} = E_{\text{POT}} + E_{\text{KIN}}$$

Its value is expected to be constant during a free fall.

The position of  $L = 0$  on the tape is important (in contrast with the acceleration measurement). Try to release the tape instantaneously.

The speed  $v$  is found like before, and the length of the fall  $L$  is measured directly on the tape:

$$v = \frac{\Delta s}{\Delta t}$$

Deviations from the conservation of energy can be assessed by calculating them as percentages of the total mechanical energy  $E_0$ .

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© Frederiksen Scientific A/S