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AC 2872.10

## This Accessory Set contains:

- 1 fluorescent red piece of plastic.
- 1 phosphorescent yellow piece of plastic.
- 1 piece of ordinary window glass  
(marked "alm. glasplade")
- 1 piece of ordinary Plastic Film  
(marked "alm. plasticfilm")
- 1 piece of UV-plastic film  
(marked "UV-plasticfilm")
- 1 UVA-filter 300-400 nm

## The different parts have the following properties:

Fluorescent red piece of Plastic:

Fluoresces in UV-light.

Phosphorescent yellow piece of Plastic:

Phosphorescent in UV-light.

Ordinary Glass:

A soda-lime glass, as is typically used in window construction. Material allows ~90% of UV and visible light to pass from 500 to 2700 nm.

## Ordinary Plastic Film:

Allows all visible and most UV light to pass.

UV-plastic film:

Absorbs UV light.

## UVA-filter:

Allows UV light with a wavelength of 300-400 nm to pass.

A good UV-lamp is strongly recommended, for instance Frederiksen's 2871.00 UV-lamp which has a band of UV-light around 370 nm (and additional 4 spectral lines in the visible range).

## Suggested Activities:

Demonstrate that the fluorescent piece of plastic still fluoresces when covered by the glass plate, the ordinary plastic film, or the UVA-filter, but loses this ability when covered by the UV-plastic film.

Demonstrate the same with the phosphorescent piece of plastic. If a UV-light is looked at, through the UVA-filter, a blue light can be seen. If a UV-light is looked at through both the UVA-filter and the UV-plastic film, no light is seen.

Apply different suntan lotions on the ordinary plastic film, and use the fluorescent piece of Plastic to observe how UV-light is absorbed by the suntan lotions.

If a more quantitative measurement is required, an UVA-sensor (like 2872.00) can be used to measure the intensity of the UVA-light.

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© A/S Søren Frederiksen, Ølgod