Vejledning til strømforsyning nr. 3610.60

06.06.12

Ac 3610.60

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06.06.12

Ac 3610.60

The power supply no. 3610.60 provides AC and DC power in 2 volt increments, from 1 to 24 volts, at up to 5 A. As the AC and DC power are delivered on separate output connectors, they may be used simultaneously.

Both output are unregulated, and protected from current overload by thermal circuit breaker. The DC output is fullwave rectified, but unfiltered.

Operation:

1. Connect the power chord to the mains outlet 230 volt AC.
2. Connect the appropriate output lines of the power supply to the circuit (test leads are not included).
3. Set the »VOLTAGE ADJUST« knob to 2 volt.
4. Flip the »On/off« switch on the power supply to On. The switch will light to show that the power supply is on.
5. Set the ?VOLTAGE ADJUST? knob to the desired voltage setting (AC values are in RMS)

Note: Both the AC and DC output circuit are protected with a thermally activated circuit breaker. If at any time the power supply fails to come on, or if it shuts down during operation because of a current overload, press the ?Reset? button on the front panel. This will reset the circuit breaker. Of course, you may have to reduce the voltage or increase the resistance of the load in order to resume operation.

Tecnical data:

Fuse: 2 x 4 A slow fuse no. 4093.13.

Total max. load: 120 watt.

Dimensions: LxWxH: 203 x 225 x 117 mm.

Weight: 2.8 kg excl. mains cable.

DC:

Voltage: 2-24 volt (unfiltered), adjustable in 2 volt increments.

Current: Max. 5 A.

Fuse: Thermal cut-out at 6 A.

AC:

Voltage: 2-24 volt (RMS), adjustable in 2 volt increments.

Current: Max. 5 A.

Fuse: Thermal cut-out.

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© A/S Søren Frederiksen, Ølgod