Strømforsyning, 0-6000 V DC nr. 3660.50

13.06.12

Ac 3660.50

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13.06.12

Ac 3660.50

The kilovolt powersupply provides regulated DC power from 0 to 6 kilovolts at up to 2 mA. The 6 kilovolt output also has a center tap, providing outputs of 0 to ± 3 kilovolts. In addition, the supply output jacks provide a constant 6.3 volts of AC power in electron tube experiments.

CAUTION - HIGH VOLTAGE:

Always work with a partner
Wear rubber soled shoes.
Work with one hand only
Keep the other hand in your pocket or behind your back.

To use the Power Supply:

1. Turn the High Voltage Adjust knob fully counterclockwise.
2. Plug the power cord into the outlet.
3. Connect the appropriate output terminals of the power supply to the circuit.
4. Flip the On/Off switch on the front panel to On. The switch will light to show that the power is on.

5. Turn the High Voltage Adjust knob to obtain the desired output. The meter will show the voltage between the 6 kV output terminals. If you are using either set of 3 kV terminals, the meter will read twice the actual output voltage.

All output circuits are protected against current overload. If at any time the power supply fails to turn on, or if it shuts down during operation because of a current overload, check the fuse on the rear of the supply. If it is blown, replace only with a similarly rated fuse.

Specifications:

Output: DC 0 to 6 kV and 0 to ± 3 kV max 2 mA.
AC - 6.3 V / 3 A

Input:

3660.50: 230 Volt AC - 50 Hz. Fuse 500 mA slow
3660.51: 115 Volt AC 50/60 Hz. Fuse 1 A slow

Warning !

To avoid electric shock the power cord protective grounding connector must be connected to ground.

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© A/S Søren Frederiksen, Ølgod