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The power meter is a versatile electric measuring instrument that apart from power also measures voltage, current, energy, charge and time.

The measured energy is summed in one of two memory banks – the same applies to the charge.

Functions

The sockets marked *In* connects to a power supply, an accumulator, a solar panel or similar.

Positive voltage goes to the red socket.

The sockets marked *Out* connects to the load like a resistor, an incandescent lamp or similar.

Voltage, current and power are all read in the usual units of V, A and W.

Energy is given in the unit Wh (watt hours). Conversion to and from joule:

$$1 \text{ Wh} = 3600 \text{ J}$$

$$1 \text{ J} = 2,7778 \cdot 10^{-4} \text{ Wh}$$

Charge is given in the unit Ah (ampere hours). Conversion to and from coulomb:

$$1 \text{ Ah} = 3600 \text{ C}$$

$$1 \text{ C} = 2,7778 \cdot 10^{-4} \text{ Ah}$$

Energy and charge are summed in memories “0” and “1”. The choice is controlled by the small pushbutton to the left of the display:

- Press button briefly to change
- Press and keep the button depressed for a few seconds to zero

Time (from turning on the instrument) is displayed in hours and minutes. The time can only be reset by turning off.

Display

- | | |
|---------------------|-------|
| 1. Time | hh:mm |
| 2. Power | W |
| 3. Memory | 0 / 1 |
| 4. Int. temperature | °C |
| 5. Charge | Ah |
| 6. Voltage | V |
| 7. Current | A |
| 8. Energy | Wh |

Current and voltage measurements

The voltage is measured across the *input sockets*. The voltage across the output is slightly lower because of the current measurement.


The current is measured in the *return path*. The instrument monitors the current running in by the black *Out* socket and out of the black *In* socket.

In most applications there is no need to pay attention to these details.

Examples of applications

Measuring energy production

This setup can for instance be used for measuring the daily energy production of a solar panel:

Connect a 488540 Solar Panel to the input of the power meter.

Adjust the tilt of the panel and make it face southwards.

The output of the power meter is connected to a 422070 Rheostat, 100 Ohm.

The rheostat must be adjusted to make the solar panel deliver maximum electric power at noon (in full sunshine). This is easy to do by observing the power reading of the meter.

If you wish to start measuring in the morning, the optimum value cannot be found experimentally. Instead, simply adjust the rheostat to 14 Ω using an ohmmeter (disconnect temporarily from the power meter). That is a suitable value.

Turn the power meter off and on to reset the time. Press the small button for three seconds to reset the energy measurement.

Wait 24 hours. Read the total energy.

Note: Be careful not to expose the power meter or the rheostat to moisture.

Characteristic curve of a solar panel

The same setup can also be used for measuring the characteristic curve of the solar panel. A halogen floodlight can substitute for the sun.

Write down voltage, current and power for a range of settings of the rheostat from short circuit (0 Ω) to maximum (100 Ω). Take care to have

readings close to the maximum power. Finally, unplug one of the leads to the rheostat and take a last reading.

Plot voltage as a function of current. This is called the solar panel's characteristic curve.

Plot also power as a function of current. Find the current where power is maximum. Mark the corresponding point on the characteristic curve.

Specifications

Voltage

Measuring range	0 – 30.00 V
Accuracy	$\pm 0.3\% + 2$ digits

Current

Measuring range	0 – 3.000 A
Accuracy	$\pm 0.8\% + 3$ digits
Voltage drop	< 200 mV (typ.)

Power

0 – 90.00 W

Energy

0 – 999.99 Wh

Time

00:00 - 99:59

(The instrument also contains a built-in, rather inaccurate thermometer.)

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© Frederiksen Scientific A/S