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4530.10 AC

Voltameter

For the demonstration of the gas quantity developed through electrolysis.

The apparatus consists of a clear plastic jar with two platinum electrodes with 4 mm terminals. Supplied with 2 testtubes to collect the developed gas.

Supplied with a movable mm scale to measure the height of the gas volume. Height 240 mm.

Electrolysis of a thin dilution of Sulphuric Acid (H_2SO_4) will generate oxygen (O_2) at the positive electrode and hydrogen (H_2) at the negative electrode. Use a thin dilution of sulphuric acid H_2SO_4 for the experiment.

Operation

1. Pour a thin dilution of sulphuric acid into the app. cup. The concentration ought not to exceed 1 M.
2. Fill the two mini test tubes with the diluted acid and place them in the app. as illustrated above.
3. Apply a voltage of 10-20 V DC to the app. terminals.

After a while gases will have accumulated in the test tubes. The gases accumulated are hydrogen (H_2) at the negative electrode and oxygen (O_2) at the positive electrode. Hydrogen and oxygen are generated in the ratio 2:1. The result of the electrolysis is a decomposition of the water in its two components hydrogen and oxygen.

Avoid electrolytes containing chlorine (e.g. NaCl solution)

Rinse the apparatus thoroughly in running water and dry it carefully before you put it aside.

Trouble shooting:

If the generation of gas proceeds slowly – try rinsing the gasvoltameter in a dilution of ammonia water. Do not use alcohol or other solvents.

NB! The gasvoltameter no. 4530.10 cannot endure exposure to ethanol even in weak dilutions.

Technical data:

Cylinder (inner dia.): 84 mm
Height (inner): 50 mm

Weight incl. test tubes: 620 g

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A/S Søren Frederiksen, Ølgod