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
The current balance consists of a magnet assembly, six interchangeable conductors and a tip-up arm with sockets for 4 mm safety plugs.

Use this equipment to study the force on a current carrying conductor in a magnetic field.

Equipment needed

Digital scales, 200 g / 0.01 g (102950, 103245, 102905 etc.). The scales must allow zeroing with the magnet assembly placed on it. The magnet assembly weighs close to 170 g. After zeroing, the scales must be able to measure both positive and negative values.

Power supply with ammeter (364000).

An adjustable current limiter is strongly recommended.

Stand (000100, 000850).

Cables (105720, 105721)

Vernier caliper or ruler.

Optional: **Teslameter** (406050) – see text.

Recommended accessory

456510 Current balance, angle-dependent

The basic current balance 456500 uses a fixed geometry. With 456510, the angle between conductor and B -field can vary more than $\pm 90^\circ$.

456510

Measurement principles

When a current I runs in a conductor placed perpendicular to a magnetic field B , the conductor experiences a magnetic force given by

$$F = B \cdot I \cdot L$$

where L is the length of the conductor. (This formula is called Laplace's force law and is a special case of Lorentz' force law.)

The magnetic field is created by 6 replaceable permanent magnets. The strength of the field is roughly (but not exactly) proportional to the number of magnets in the holder.

The precise value can be measured with e.g. our teslameter 406050. (In many cases it is sufficient to know that B is a constant.)

The conductor is a track on a printed circuit board that is mounted at an arm that can be tilted for easy exchange of circuit boards.

There are six different circuit boards with different conductor lengths. This is measured as the length of the *horizontal* part of the track; The two vertical sections are acted on by equally large forces with *opposite* sign.

The current is measured with an ammeter (possibly included in the power supply). It is recommended to use a power supply with a variable current limiter as the low resistance in the circuit makes it difficult to keep the current at a given value.

If you don't have this kind of power supply, a 1Ω 50 W resistor must be inserted in series with the current balance (use e.g. 429140).

The force on the magnet has the same magnitude as the force on the conductor, but the opposite direction (Newton's 3rd law: Action = reaction). The magnet assembly is placed on a scale, the reading of which can be used to find the force on the conductor:

$$F = m \cdot g$$

where m is read on the scale and g is the acceleration due to gravity.

Experiments

Frederiksen Scientific has produced the following complete lab manuals, available for free at www.frederiksen.eu:

137230 The magnetic force on a conductor
(uses 456500)

137240 Laplace's force law
(uses 456500 and 456510)

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© Frederiksen Scientific A/S