2018-05-29 / HS

Produktmanual AC 489030

English
manual
overleaf

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Product manual AC 489030

Dansk
manual:
Vend!

The student pyranometer is an affordable instrument for measuring solar radiation on a surface.

The instrument connects to an ordinary voltmeter (multimeter).

The pyranometer is temperature compensated.

Using the pyranometer

Connect the pyranometer sockets to a multimeter with a pair of lab leads. Select DC voltage. The measurement range should be chosen to allow precise readings around 100 mV. The pyranometer is calibrated to output a voltage of 100 mV when placed perpendicular to solar radiation with an intensity of 1000 W/m². This corresponds approximately to the intensity of sunlight on a cloudless day.

To find the radiant flux on a given surface (e.g. a solar collector), place the pyranometer parallel to the surface. Read the voltmeter and convert to power per area. Knowing the area of the surface enables you to calculate the total radiant flux.

Example

A small solar collector model is placed perpendicular to the sunbeams. The student pyranometer is placed facing the same direction. The voltmeter reading is 94.5 mV. That means that the radiant flux density is 945 W/m².

If the solar collector measures 25 x 25 cm (i.e. 0.0625 m²), the total flux (or power) amounts to

$$945 \text{ W/m}^2 \cdot 0.0625 \text{ m}^2 = 59 \text{ W}$$

Specifications

These pyranometers are neither calibrated nor controlled individually. The figures below are *typical* data based on sample inquiries.

Sensitivity: 0.10 V / (W/m²)

Absolute calibration: ±5 %

Temperature variations: <±2 % (-10°C to +75°C)

The specified sensitivity applies to sunlight. Using light sources with another spectral distribution will result in a slightly different sensitivity.

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© Frederiksen Scientific A/S