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# Manual for Electric Motor with Propeller no. 5015.00

20.12.10

AC 5015.00

This motor distinguish itself by its extremely low energy consumption, which makes it especially well suited for

demonstration of low voltage/current equipment e.g. thermal energy produced by the Thermogenerator no.

4875.00 or photovoltaic energy produced by a Solar Cell as e.g. no. 4885.00. The propeller will start rotating at

as low an input voltage as 35 mV! Please refer to the technical data below.


N.B. Please notice that the propeller axis is very thin and delicate and thus susceptible to mechanical action.

## Motor data

1. Nominal voltage . . . . .  $U_N$  . . . . . 1,5 Volt
2. Recommended speed . . . . .  $N_E$  . .  $\leq 12.000$  rpm
3. Recommended torque . . . . .  $M_E$  . . . .  $\leq 0,5$  mNm
4. Recommended power . . . . .  $P_{2E}$  . . . . .  $\leq 0,3$  W
5. Recommended current . . . . .  $I_E$  . . . . .  $\leq 0,596$  A

## Lemon energy

Quite a few illustrations concerning production of chemical energy shows a lemon with a zinc and a copper electrode inserted. This simple cell is connected to an electric motor with propeller as e.g. no. 5015.00. In reality the propeller doesn't always start to rotate once the electrodes has been inserted into the lemon. If you study the matter more closely it is possible to make the arrangement work. You have to optimise things to an extent that doesn't always quite match the result. The crucial points are

a careful cleaning of the electrodes by means of steel wool and/or ammonia water, the distance between the electrodes and the condition of the lemon - the fresher the better. A simpler and better way to demonstrate these matters are as follows: Soak a piece

of filter-paper approx. 10 cm dia. with lemon juice from a fresh lemon, place two clean electrodes on each side of the filter-paper as illustrated. This system is in fact a simple cell. Connect the electrodes to the Electric Motor no. 5015.00 and the propeller will start rotating. Should it fail to rotate or if it stops after a while then simply remove one of the electrodes from the filter-paper and give it an airing for approx. 5-10 sec. Press the electrode against the filter-paper again and the propeller will start rotating. The propeller will stop rotating again after a while as the chemical process will produce Hydrogen which obstructs the process. Remove one of the electrodes again and allow the hydrogen to vanish. The procedure may be repeated over and over again. This method also offers you an excellent option to delve into the chemical process. The electrodes should be cleaned carefully by means of steel wool and subsequently degreased by means of orange e.g. ammonia water.

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© A/S Søren Frederiksen, Ølgod