QUICK-GUIDE

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QUICK GUIDE

Power on: Press ON to power on the instrument. The display briefly shows a version number and the GM voltage. Immediately after the counter has been turned on, the counting period is 1 s – **CPS** is indicated on the display.

Other views: From CPS press ▲ once to select **µSv/h** relative to Cs-137.

Press ▲ once more to select the mean value of counts per second – **Mean CPS**.

Other counting periods: Press ▲ or ▼ until **Pulse** is highlighted in the display. The instrument now counts for a selectable period, starting with 10 s.

Press HOLD ▲ HOLD and set the measurement period with ▲ and ▼.

The unit is seconds. The symbol “---” means indefinite counting.

Select with **CAL** / cancel with **HOLD**.

Setting the GM voltage: Press **CAL** and adjust the voltage with ▲ and ▼. Select with **CAL**.

Counting rate (**CPS**) is shown in the secondary display during the adjustment.

This setting returns to the default value when the instrument is turned on.

- 8 fixed count periods (1, 10, 20, 30, 40, 50, 60 or 100 s)
- Display of dose rate ($\mu\text{Sv}/\text{h}$) related to Cs-137
- GM tube voltage can be adjusted between 400 V and 700 V while showing count rate.
- Max, Min and average view for $\mu\text{Sv}/\text{h}$ and CPS.

Description

The GM tube

The Geiger tube is provided with a cap to protect the sensible window. When working with gamma radiation, the cap may remain in place. Remove the cap when alpha or beta radiation is involved.

Time

The counting period 1 s is achieved by selecting **CPS** (counts per second). All other periods, including infinity is obtained with the **Pulse** setting. When the instrument has been turned off, the measurement period for **Pulse** is 10 s.

To change the period, press **HOLD ▲ HOLD** and set the measurement period with ▲ and ▼.

The unit is seconds. The symbol “---” means infinite counting.

Select with **CAL** / cancel with **HOLD**.

During Pulse measurements, the progress of each measurement can be followed on the % scale at the top

GM voltage

The tube voltage can be adjusted in increments of 10 V. Nominal voltage is 500 V.

It is possible to follow the count rate (**CPS**) in the secondary display while the voltage is adjusted. Thus one can fairly easily estimate the knee voltage of the tube. This should be done by *increasing* the voltage as the response time is somewhat longer when the voltage is reduced.

Press **CAL** and adjust the voltage with **▲** and **▼**. Select with **CAL**.

This setting returns to default value when the instrument is turned off.

Setting the default value for the GM voltage

Switch off with OFF. Keep **CAL** pressed while turning on again (until the symbol **CAL** is displayed).

Set the voltage with **▲** and **▼**. Select with **CAL**.

Statistics

The instrument keeps track of the maximum, minimum and average value of the **CPS** and **µSv/h**. The statistics functions starts when the instrument is switched on.

The values are displayed like this:

| | | | | |
|-----------------|----------------|------------|--|--------------------|
| HOLD ▼ | Display: [HOLD | MAX | shows max. CPS and µSv/h | HOLD = exit |
| HOLD ▼▼ | Display: [HOLD | MIN |] Shows min. CPS and µSv/h | HOLD = exit |
| HOLD ▼▼▼ | Display: [HOLD | AVG |] Shows avg. CPS and µSv/h | HOLD = exit |

Resetting

Instead of pressing **HOLD** in the end, press **CAL** to restart the statistics functions and use the current values for **CPS** and **µSv/h** as new starting point.

The instrument has a USB-like socket placed under a cover. This is *not* for connection of USB devices.
(It is used during production for test purposes.)

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© Frederiksen Scientific A/S