The Biotechnology Education Company ©

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[www.edvotek.com](http://www.edvotek.com)

The PCR process and Taq DNA polymerase are covered by patents owned by Hoffman-LaRoche, Inc.

345.130821

EDVOTEK, The Biotechnology Education Company, and InstaStain are registered trademarks of EDVOTEK, Inc. UltraSpec-Agarose, PCR EdvoBead, and FlashBlue are trademarks of EDVOTEK, Inc.

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Experiment # 345 contains material for up to 25 reactions.

Sample volumes are very small. It is important to quick spin the tube contents in a microcentrifuge to obtain sufficient volume for pipetting. Spin samples for 10-20 seconds at maximum speed.

All components are intended for educational research only. They are not to be used for diagnostic or drug purposes, nor administered to or consumed by humans or animals.

### Materialer

	Materialer	Opbevaring	Check (✓)
A	Tubes with PCR EdvoBeads™ Each PCR EdvoBead™ contains <ul style="list-style-type: none"> <li>• dNTP Mixture</li> <li>• Taq DNA Polymerase Buffer</li> <li>• Taq DNA Polymerase</li> <li>• MgCl<sub>2</sub></li> <li>• Reaction Buffer</li> </ul>	Room Temperature	<input type="checkbox"/>
B	PTC Primer mix concentrate	-20°C Freezer	<input type="checkbox"/>
C	100 base pair ladder	-20°C Freezer	<input type="checkbox"/>
D	Control DNA concentrate	-20°C Freezer	<input type="checkbox"/>
E	TE buffer	-20°C Freezer	<input type="checkbox"/>
F	Proteinase K	Room temperature	<input type="checkbox"/>
G	Restriction Enzyme Reaction Buffer	-20°C Freezer	<input type="checkbox"/>
H	Restriction Enzyme Dilution Buffer	-20°C Freezer	<input type="checkbox"/>
I	Hae III Restriction Enzyme <ul style="list-style-type: none"> <li>• PTC Paper</li> <li>• Control Taste Paper</li> </ul>	Room temperature Room temperature	<input type="checkbox"/> <input type="checkbox"/>

**NOTE:** Components B and D are supplied in concentrate form and require dilution prior to setting up PCR reactions.

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
### OPTIONAL STOPPING POINT:

The extracted DNA may be stored at -20°C for amplification at a later time.

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
- PLACER** rørene i is

**PROCEED** to Module III: Restriction Digest of PCR Product.

**OPTIONAL STOPPING POINT:** The PCR samples may be stored at -20° C for restriction digest at a later time.

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
**PROCEED** to Module IV: Agarose Gel Electrophoresis.

**OPTIONAL STOPPING POINT:** The restriction digests may be stored at -20°C for electrophoresis at a later time.

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

Table

A

### 1x Electrophoresis Buffer (Chamber Buffer)



EDVOTEK Model #	Total Volume Required	Dilution	
		50x Conc. Buffer	+ Distilled Water
M6+	300 ml	6 ml	294 ml
M12	400 ml	8 ml	392 ml
M36	1000 ml	20 ml	980 ml

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Table  
B

Individual 2.0%  
UltraSpec-Agarose™ Gel

Size of Gel Casting tray	1X Diluted Buffer	+ Amt of Agarose
7 x 7 cm	25 ml	0.50 g
7 x 14 cm	50 ml	1.0 g

Table  
C

Time and Voltage Guidelines  
(2.0% Agarose Gels)

Volts	Time: 7 x 7 cm gel ~4.0 cm migration	Time: 7 x 14 cm gel ~6.5 cm migration
125	30 min.	60 min.
70	60 min.	120 min.
50	90 min.	150 min.

Table 1

Gel 1: PCR Reaction for PTC Gene	
Lane 1	100 bp ladder
2	Undigested PCR Product - Control
3	Undigested PCR Product - Taster 1
4	Undigested PCR Product - Taster 2
5	Undigested PCR Product - Taster 3
6	Undigested PCR Product - Taster 4
Gel 2: RFLP Analysis to Detect PTC Polymorphisms	
Lane 1	100 bp ladder
2	<i>HaeIII</i> Digested PCR Product - Control
3	<i>HaeIII</i> Digested PCR Product - Taster 1
4	<i>HaeIII</i> Digested PCR Product - Taster 2
5	<i>HaeIII</i> Digested PCR Product - Taster 3
6	<i>HaeIII</i> Digested PCR Product - Taster 4

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#### FORBEREDELSE INDEN DE ENKELTE MODULER:

Preparation For:	What to do:	When:	Time Required:
<b>Module I: Isolation of DNA from Cheek Cells</b>	Prepare and aliquot various reagents (Saline, Lysis buffer)	Up to one day before performing the experiment. <b>IMPORTANT:</b> Prepare the Lysis buffer no more than one hour before performing the experiment.	30 min.
	Equilibrate waterbaths at 55° C and boiling.	One hour before performing the experiment.	5 min.
<b>Module II: Amplification of the PTC Regions</b>	Prepare and aliquot various reagents (Primer, DNA template, ladder, etc.)	One day to 30 min. before performing the experiment.	30 min.
	Program Thermal Cycler	One hour before performing the experiment.	15 min.
<b>Module III: Restriction Digest of the PTC PCR Products</b>	Prepare and aliquot various reagents (Reaction buffer, <i>HaeIII</i> Restriction enzyme, etc.)	One day to 30 min. before performing the experiment. <b>IMPORTANT:</b> Prepare the diluted <i>HaeIII</i> Restriction Enzyme no more than one hour before performing the experiment.	30 min.
	Equilibrate waterbath at 37° C	One hour before performing the experiment.	15 min.
<b>Module IV: Separation of PCR Products &amp; Digestion Products by Electrophoresis</b>	Prepare diluted TAE buffer	Up to one day before performing the experiment.	45 min.
	Prepare molten agarose and pour gel		
<b>Module V: Staining Agarose Gels</b>	Prepare staining components	The class period or overnight after the class period.	10 min.
<b>Module VI: Confirmation of Bitter Tasting Ability</b>	Distribute PTC taste strips and Control taste strips	The class period	5 min.

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## Appendices

- A EDVOTEK® Troubleshooting Guide
- B Preparation and Handling of PCR Samples With Wax
- C Bulk Preparation of Agarose Gels

Material Safety Data Sheets:

Now available for your convenient download on [www.edvotek.com](http://www.edvotek.com).

EDVO-TECH Service

**1.800.EDVOTEK**

(1.800.330.6035)

**Please Have the Following Info:**

- Experiment number and title
- Kit lot number on box or tube
- Literature version (in lower right corner)
- Approx. purchase date

Mon-Fri  
8am-5:30pm ET

FAX 202.370.1501 • [info@edvotek.com](mailto:info@edvotek.com) • [www.edvotek.com](http://www.edvotek.com)

Visit Us Online!

**[www.edvotek.com](http://www.edvotek.com)**

**Download**  
Edvotek kits  
with  
Protocols!

**Browse**  
our online  
catalog!

**Order**  
**Products**  
online!

**Receive**  
Technical  
Support!

EDVOTEK - The Biotechnology Education Company®  
1-800-EDVOTEK • [www.edvotek.com](http://www.edvotek.com)  
FAX: 202-370-1501 • email: [info@edvotek.com](mailto:info@edvotek.com)

345.130821

## EDVOTEK® Troubleshooting Guides

## DNA EXTRACTION

PROBLEM:	CAUSE:	ANSWER:
There is no cell pellet after centrifuging the cheek cell suspension.	Not enough cheek cells in suspension	Mouth must be vigorously rinsed for at least 60 sec to harvest loose cheek cells.
	Sample not centrifuged fast enough	Spin cells at maximum speed (17,000 x g) for 2 min. If your centrifuge does not reach this speed, spin at highest available speed for 4 min.
Poor DNA Extraction	Samples not mixed well enough during extraction	In addition to flicking the tube, vortex or pipet up and down to mix the sample.
	Proteinase K inactive because it was prepared too far in advance.	Prepare Proteinase K within one hour of use.
	Water baths not at proper temperature	Use a thermometer to confirm water bath set point.
	Not enough DNA	Try cheek cell extraction. Final DNA concentrations are usually higher.
The extracted DNA is very cloudy.	Cellular debris from pellet transferred to tube	Centrifuge sample again and move supernatant to a fresh tube. Take care to avoid pellet.
	Cellular debris not separated from supernatant	Centrifuge sample again. If possible, centrifuge at a higher speed. Move cleared supernatant to a fresh tube.

## RESTRICTION ENZYME DIGESTION

PROBLEM:	CAUSE:	ANSWER:
Undigested or incompletely digested DNA	Impure DNA – some contaminants (EDTA, salts) might partially or completely inhibit activity of <i>Hae</i> III restriction enzyme	Poor DNA extraction. Extract new DNA. Cheek cell extraction usually results in higher DNA yield.
	Improper dilution of enzyme	Ensure that <i>Hae</i> III restriction enzyme was correctly diluted.
	Improper addition of enzyme	Ensure that correct amount of <i>Hae</i> III restriction enzyme was added to the restriction digest.
	Incorrect incubation temperature	Use a thermometer to confirm water bath temperature and adjust, if necessary.
Unexpected Cleavage Pattern	DNA sample is contaminated	Prepare a new DNA sample.
Smearing of digested DNA on gel	Nuclease contamination	Care should be taken to avoid cross contamination when setting up reactions.
	Agarose running conditions	Use fresh electrophoresis buffer and appropriate voltage.

## EDVOTEK® Troubleshooting Guides

## PCR AND ELECTROPHORESIS

PROBLEM:	CAUSE:	ANSWER:
There is very little liquid left in tube after PCR	Sample has evaporated	<p>Make sure the heated lid reaches the appropriate temperature.</p> <p>If your thermal cycler does not have a heated lid, overlay the PCR reaction with wax (see Appendix B for details)</p> <p>Make sure students close the lid of the PCR tube properly.</p>
	Pipetting error	Make sure students pipet 20 $\mu$ L primer mix and 5 $\mu$ L extracted DNA into the 0.2 mL tube.
The ladder, control DNA, and student PCR products are not visible on the gel.	The gel was not prepared properly.	<p>Ensure that the electrophoresis buffer was correctly diluted.</p> <p>Gels of higher concentration (&gt; 0.8%) require special attention when melting the agarose. Make sure that the solution is completely clear of "clumps" and glassy granules before pouring gels.</p>
	The gel was not stained properly.	Repeat staining.
	Malfunctioning electrophoresis unit or power source.	Contact the manufacturer of the electrophoresis unit or power source.
After staining the gel, the DNA bands are faint.	The gel was not stained for a sufficient period of time.	Repeat staining protocol.
After staining, the ladder is visible but no PCR products are present.	PCR amplification was unsuccessful.	Repeat PCR with fresh PCR EdvoBeads™ and primers.
		Ensure that the thermal cycler has been properly programmed. See Module II for guidelines
After staining, the ladder and control PCR products are visible on gel, but some student samples are not present.	Student DNA sample was not concentrated enough.	Poor DNA extraction. Extract new DNA. Cheek cell extraction usually results in higher DNA yield.
	Student DNA sample was degraded	If DNA is not used immediately following extraction, store sample at -20°C.
	Wrong volumes of DNA and primer added to PCR reaction	Practice using pipettes
Some students have more or less amplification than others.	Concentration of DNA varies by sample.	There is an inherent variability in the extraction process. For best results, use cheek cell extraction.
Low molecular weight band in PCR samples	Primer dimer	Low concentration of extracted DNA in PCR reaction.
DNA bands were not resolved.	Tracking dye should migrate at least 3.5 cm (if using a 7x7 cm tray), and at least 6 cm (if using a 7x14 cm tray) from the wells to ensure adequate separation.	Be sure to run the gel at least 6 cm before staining and visualizing the DNA (approximately one hour at 125 V).
DNA bands fade when gels are kept at 4°C.	DNA stained with FlashBlue™ may fade with time	Re-stain the gel with FlashBlue™

### Preparation and Handling of PCR Samples With Wax

ONLY For Thermal Cyclers WITHOUT Heated Lids, or Manual PCR Using Three Waterbaths

Using a wax overlay on reaction components prevents evaporation during the PCR process.

#### HOW TO PREPARE A WAX OVERLAY

1. Add PCR components to the 0.2 ml PCR Tube as outlined in Module III.
2. Centrifuge at full speed for five seconds to collect sample at bottom of the tube.
3. Using clean forceps, add one wax bead to the PCR tube.
4. Place samples in PCR machine and proceed with Module III.

#### PREPARING PCR SAMPLES FOR ELECTROPHORESIS

1. After PCR is completed, melt the wax overlay by heating the sample at 94° C for three minutes or until the wax melts.
2. Using a clean pipet, remove as much overlay wax as possible.
3. Allow the remaining wax to solidify.
4. Use a pipet tip to puncture the thin layer of remaining wax. Using a fresh pipet tip, remove the PCR product and transfer to a new tube.
5. Add 5  $\mu$ L of 10x Gel Loading Buffer to the sample. Proceed to Module IV to perform electrophoresis.

## Bulk Preparation of Agarose Gels

To save time, electrophoresis buffer and agarose gel solution can be prepared in larger quantities for sharing by the class. Unused diluted buffer can be used at a later time and solidified agarose gel can be remelted.

## BULK ELECTROPHORESIS BUFFER

Bulk preparation of 1X electrophoresis buffer is outlined in Table D.

Table  
D

Bulk Preparation of Electrophoresis Buffer

50x Conc. Buffer	+	Distilled Water	Total Volume Required
60 ml		2,940 ml	3000 ml (3 L)


## BATCH AGAROSE GELS (2.0%)

Bulk preparation of 2.0% agarose gel is outlined in Table E.

1. Use a 500 ml flask to prepare the diluted gel buffer.
2. Pour the appropriate amount of UltraSpec-Agarose™ into the prepared buffer. Swirl to disperse clumps.
3. With a marking pen, indicate the level of solution volume on the outside of the flask.
4. Heat the agarose solution as outlined previously for individual gel preparation. The heating time will require adjustment due to the larger total volume of gel buffer solution.
5. Cool the agarose solution to 60°C with swirling to promote even dissipation of heat. If evaporation has occurred, add distilled water to bring the solution up to the original volume as marked on the flask in step 3.
6. Dispense the required volume of cooled agarose solution for casting each gel. The volume required is dependent upon the size of the gel bed.
7. Allow the gel to completely solidify. It will become firm and cool to the touch after approximately 20 minutes. Proceed with electrophoresis (Module IV) or store the gels at 4°C under buffer.

**Note:**

The UltraSpec-Agarose™ kit component is usually labeled with the amount it contains. Please read the label carefully. If the amount of agarose is not specified or if the bottle's plastic seal has been broken, weigh the agarose to ensure you are using the correct amount.

Table  
EBatch Prep of 2.0%  
UltraSpec-Agarose™

Amt of Agarose	+	Diluted Buffer (1x)
8.0 g		400 ml

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© 2013 Edvotek® All Rights Reserved.